

MINI MOBILE WITH FALL DETECTION

USER GUIDE

HOW TO CHARGE YOUR MINI MOBILE WITH FALL DETECTION

IMPORTANT:

Fully charge your Mini Mobile with Fall Detection before first use and when the red battery light is flashing while off the cradle to ensure the device can function properly.



Plug in the charging cradle and place the Mini Mobile with Fall Detection into the cradle as shown. The Mini Mobile with Fall Detection will state "Charging". Typical charging time is 2-3 hours. The red battery light flashes every five seconds while charging and is solid red when the Mini Mobile with Fall Detection is fully charged.

NOTE:

Only charge as needed. Wear your device as often as possible to be protected. If the battery is low, the device will state, "Battery low, please charge" or "Call failed, battery low" when you press the call button. It will then power off.



MEET YOUR MINI MOBILE WITH FALL DETECTION



- | | | |
|--|----------------|------------------|
| 1. Call Button | 3. Speaker | 5. Charging Pins |
| 2. Blue Call Light / Red Battery Light | 4. Microphones | |

RED BATTERY LIGHT AROUND CALL BUTTON

- While charging, the red battery light will flash every five seconds until the device is fully charged.
- Once fully charged, the red battery light will be solid.
- If the battery is low, the red battery light will slowly flash when off the charging cradle.

BLUE CONNECTION LIGHT AROUND CALL BUTTON

- During an SOS call, the blue light will be solid.
- The blue light will also flash every five seconds while charging.
- The light may flash blue every 26 seconds when off the cradle, depending on operating mode.
- If the blue light flashes every second, or if it double flashes, please contact Galaxy Medical Alert at 1-855-905-7414.

HOW TO CALL FOR HELP WITH THE MINI MOBILE WITH FALL DETECTION



1. Press and hold the call button and count 1, 2,3, or until the light around the call button on the Mini Mobile with Fall Detection turns blue, then release the button.
2. After a short delay, you will hear a voice message and then tones or ringing.
3. The specialist will answer the call.
4. The specialist can send emergency personnel or loved ones to help you.

* If your Mini Mobile with Fall Detection does not connect properly, please contact Galaxy Medical Alert at 1-855-905-7414 immediately.

HOW TO TURN OFF THE MINI MOBILE WITH FALL DETECTION

1. Turn over the Mini Mobile with Fall Detection device and locate the recessed black button above or below the charging pins.
2. Depress the black button for 3 - 5 seconds until the device informs you that it is powering OFF.

We recommend that you test the Mini Mobile with Fall Detection weekly.

1. To test, press and hold the call button and count 1, 2, 3, or until the light turns blue, then release the button.
2. After a short delay, you will hear a voice message and then tones or ringing.
3. The specialist will answer the call.
4. Tell the specialist you are testing.

* If your Mini Mobile with Fall Detection does not test properly, please contact Galaxy Medical Alert at 1-855-905-7414 immediately.

REPLACE THE WRISTBAND

1. Turn the Mini Mobile with Fall Detection over and locate the small silver quick release pins on each wrist strap.
2. To remove the existing wrist straps, slide the quick release pin toward the center of the device and pivot the strap away from the device.
3. To put the new wrist strap on the device, insert one side of the horizontal bar on the new band into the device, and then slide the quick release pin toward the center of the device.
4. Pivot the strap into the device, sliding the quick release pin back into position once the strap is firmly in the device.



CLEANING THE MINI MOBILE WITH FALL DETECTION WEEKLY

We recommend that you clean your Mini Mobile with Fall Detection weekly to ensure proper charging. Take a soft cloth, such as an eyeglass cleaning cloth, and gently rub the gold contacts on the device and charging cradle to remove dirt, oil and debris. You can also spray household cleaner onto the cleaning cloth (though do not spray directly onto the device or cradle) to remove debris.

Example of a customer cleaning the gold charging contacts on the bottom of the Mini Mobile with Fall Detection.



Example of a customer cleaning the speaker holes on the top of the Mini Mobile with Fall Detection.

